

Nonwoven Plaster

Technical data sheet

**ID name:** Non woven

Product qualities: Good air permeability.

Intended use:

For use on minor cuts and abrasions.
Protects minor wounds.

Direction for use:

Clean wound and surrounding area and dry skin before applying the plaster.
Position pad over wound and press firmly on adhesive area to ensure adhesion.
Change dressing every 24 hours or if it is wet.
Keep dry and out of direct sunlight.
Store below 25°C.

Caution for use:

Do not use with major wounds.
If wound bleeds excessively or skin becomes irritated or infected remove plaster and consult a doctor immediately.
If you are diabetic consult your doctor before use.

CE Class: I**Product properties:****Absorbent pad**

Raw Material: 70% Viscose+30% Polyester
Weight per unit area: 110gsm
Absorptive capacity ≥ 600%

Film material:

Material: Nonwoven
Size: 10x6cm

Adhesive coating

Henkel Acrylate adhesive

Glassing Paper

Weight per unit area: 60 gsm

Bioburden Control

Total bacteria < 100 CFU/g
Yeast / Mustiness ≤ 20 CFU/g
Pathogenic: 0

Storage:

Should be stored in a cool dry environment for best results of product performance. Keep away from light and heat (store below 25°C)

Shelf life:

5 Years

Legal manufacturer:

RFX+CARE Manufacturing Co., Ltd. 7 Lanjiang Road, Yuecheng District, Shaoxing, Zhejiang, China 312000.

EC-Representative:

RFX+Care International A/S. Bakkegaardsvej 408, 3050 Humlebaek, Denmark

Certificates TÜV SÜD:

EN ISO 13485:2016

Reference reports and certificates:

Annex. 1 ISO13485:2016 Certificate
Annex. 2 Skin Irritation Test
Annex. 3 Skin Sensitization Test
Annex. 4 Cytotoxicity Test
Annex. 5 Microbial Limit Test
Annex. 6 Conforms to DIN 13164:2014

For further information please send your request to info@rfx-care.com

EC-Representative: RFX+Care International A/S. • Bakkegaardsvej 408 • 3050 Humlebaek • Denmark