

# KLORRENT

Page 1/1 Data Sheet no.: 1055

Date: 2. jul 2012.

<b>PRODUCTTYPE</b>	Disinfectant								
<b>APPLICATION</b>	KLORRENT is a disinfectant for all chlorine resistant surfaces such as instruments, tools, machines, equipment and the like. KLORRENT is also well suited for bleaching of dish and floor cloths.								
<b>PROPERTIES</b>	KLORRENT is a bleaching agent and disinfectant. The product is used after previous cleaning.								
<b>CHARACTERISTIC</b>	Clear, light yellow liquid.								
<b>DIRECTIONS FOR USE AND DOSING</b>	For disinfection use 1-2 dl (4-8 pump strokes) for 8 liters of water. After final disinfection, all surfaces that will later be in contact with foods, shall be rinsed thoroughly with drinking water. For bleaching use 1/2 dl (2 pump strokes) for 8 liters of water. DO NOT MIX WITH ACID OR ACID CONTAINING PRODUCTS.								
<b>PRODUCT DATA</b>	<table><tr><td>Solubility</td><td>Good - dissolved immediately in water.</td></tr><tr><td>Specific volume</td><td>1,10 kg/l.</td></tr><tr><td>pH-value (diluted) app.:</td><td>11 (1%)</td></tr><tr><td>pH-value (conc.) app.:</td><td>12,5</td></tr></table>	Solubility	Good - dissolved immediately in water.	Specific volume	1,10 kg/l.	pH-value (diluted) app.:	11 (1%)	pH-value (conc.) app.:	12,5
Solubility	Good - dissolved immediately in water.								
Specific volume	1,10 kg/l.								
pH-value (diluted) app.:	11 (1%)								
pH-value (conc.) app.:	12,5								
<b>STORAGE</b>	Keep in a cool place in closed original container.								
<b>APPROVAL</b>	The disinfectant is approved by the Danish Veterinary and Food Administration under file number: 2007-20-5409-00997.								
<b>SAFETY CONDITIONS</b>	Refer to safety data sheet.								

## NOVADAN

Platinvej 21, DK-6000 Kolding, Tlf. +45 76 348 400  
PL - ul. Holenderska 4, 05-152 Czosnów, Tel +48 (0)22 785 1002